


# SEQURE<sup>®</sup>

Reflux Control Microcatheter

## Designed to Help Reduce the Reflux of Embolic Material<sup>1</sup>

The preclinical study showed that the SeQure<sup>®</sup> reflux control microcatheter delivered more microsphere volume while reducing the risk of non-target embolization.<sup>1</sup>




Guerbet |

\*Preclinical study was sponsored by Guerbet Group.  
Study authors: Maleux, De Baere, Sun, Holtzman, Harbater, Balvay, Veil, Rizzitelli, Violas, Robert, Robic, Seron, Corot.

## Fluid Barrier Technology for Flow-Directed Embolization

The SeQure<sup>®</sup> reflux control microcatheter with distal micro side slits utilizes flow dynamic principals to create a filtered contrast media fluid barrier and direct beads to the target vessel, allowing interventional radiologists to have more control over their embolization procedures.



### Reduced Non-Target Embolization and Increased Targeted Delivery with SeQure<sup>®</sup> Microcatheter Demonstrated in Swine Model Study<sup>1</sup>

#### STUDY OBJECTIVES

- A** To evaluate the performance of the SeQure<sup>®</sup> reflux control microcatheter delivering particles to the target area versus non-target area in vitro and in vivo on a swine pre-clinical model.
- B** To visualize and quantify the amount of microspheres delivered in the target site versus the non-target site; calculate the target/non target ratio.
- C** To compare the performances of SeQure<sup>®</sup> reflux control microcatheter versus a standard end-hole catheter.

#### STUDY LIMITATIONS

The use of an automatic injector for the microspheres' administration allowed the standardization of the injection but is not representative of the clinical practice. A study with a manual injection closer to clinical experience will follow this first preclinical test.

Results of the preclinical study were based on a swine model which presents a limitation in the translational aspect of the conclusion. Further exploration in a pathological model is needed.

#### STUDY METHODS

The target vessel in the kidneys of thirteen pigs were identified using 2.4 Fr (n = 7) or 2.7 Fr (n = 6) microcatheter (SeQure<sup>®</sup> versus Standard End-Hole Catheter). Evans Blue solution was injected to stain the selective area of injection followed by the injection of 75/150 µm radiopaque microspheres suspended in 50/50 lobitridol 300/NaCl solution (1 mL of beads in 10 mL of solution). Injection was performed by an automatic injector with flow rate of 8 mL/min. End point was obtained via reflux of contrast in non-target vessel at the second visible time (according to clinician experience). Angiography was performed with a 4/5 Fr catheter of the whole kidney post embolization.

Ex-vivo analysis was completed using µCT ex vivo for quantitative analysis of target embolization, non-target embolization, ratio target/non-target embolization; Evans Blue detection was for determining the target area of injection.

For more information, please contact us at 1-877-729-6679 or [medical.liason.na@guerbet.com](mailto:medical.liason.na@guerbet.com)

SeQure® is registered trademark of Guerbet Group or its affiliates.


**CAUTION:** Federal law restricts this device to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

**INDICATIONS FOR USE:** The SeQure® microcatheters are intended for the infusion of contrast media into all peripheral vessels. The SeQure® microcatheters are also intended for drug infusion in intra-arterial therapy and infusion of embolic materials. The SeQure® microcatheters should not be used in cerebral vessels.

1. Maleux et al. Reduced Non-Target Embolization and Increased Targeted Delivery with SeQure® Microcatheter Demonstrated in Swine Model. Presented at RSNA 2019.
2. Maleux et al. Reduced Non-Target Embolization & Better Delivery with an Innovative Reflux Control Microcatheter: Mechanisms of Action and Clinical Application in Four Swine. Presented at CIRSE 2019.


## Preliminary Data from the Preclinical Study<sup>2</sup>

Preliminary results on the right demonstrate that the SeQure® reflux control microcatheter showed higher node filling when compared to a standard end-hole microcatheter.



## Preclinical Results from the Preclinical Study<sup>1</sup>

Preclinical results below demonstrate that the SeQure® reflux control microcatheter produced a higher volume of microsphere (MS) delivery in the target area (and increased microspheres injected) when compared to a standard end-hole microcatheter.



## STUDY OUTCOMES

The photos below demonstrate that the SeQure® reflux control microcatheter produced better **filling of the main vessels, higher volume of microspheres injected** (green circle), and **less non-target embolization** (red box). **Ex-vivo marking** (orange circle) demonstrates more targeted fill area with the SeQure® reflux control microcatheter versus a standard end-hole microcatheter.

- Better filling of the main vessels
- Higher volume of microspheres injected
- Less non-target embolization
- More distal penetration of microspheres
- More efficient microsphere delivery
- Predictability of the treatment

