

OptiVantage[®] Single-use
power injector & its
certified* consumables
work better together

Choose the ideal
COMBINATION

* Certified by Guerbet



Essential certified by Guerbet

Combining your **OptiVantage®** with **Guerbet consumables ensures optimum performance.**

Our priority is to support you in your daily complex challenges as well as in your constant need to improve quality and to optimize your budget.

To swiftly meet your needs and to adapt to rapid market developments, we have developed trusted partnerships with selected manufacturers to provide a comprehensive portfolio of essential consumables. At all steps, we ensure that high quality standards are followed in order to provide for safe and efficient injections.



All consumables are selected, tested and certified to meet the highest Guerbet quality standards. Legal manufacturers are ISO 13485** accredited and all consumables are CE marked and comply with the European Union regulatory requirements for Medical Devices*** or with the U.S. regulations or with local regulatory requirements when products are approved or registered.

* Certified by Guerbet as compatible on Guerbet injectors - data on file. ** Legal manufacturers Liebel-Flarsheim Company LLC and Shenzhen Boon Medical Supply Co., Ltd. are ISO 13485 certified. *** The EU Medical Device Directive (93/42/EEC, MDD) is a legislation which addresses essential principles for the safety and performance of medical devices during their lifecycle. The new EU Medical Devices Regulation ((EU) 2017/745, MDR) will repeal the existing MDD and bring significant regulatory changes.



Syringe Kit

Dual syringe kit with coiled Y-tubing with dual check valve and spikes

Single-use – 200mL – 150cm – pressure limit 350psi

REF #236302



Syringe Kit for use with Prefilled Syringe

Single syringe kit with coiled Y-tubing (dual check valve) and spike


Single-use – 200mL – 150cm – pressure limit 350psi

REF #236300



Reference number	Designation article	Quantity in box	Requested quantity	Comments
236302	Dual CT Syringe Kit	20		
236300	Single CT Syringe Kit	50		
223735	150 cm Y-tubing (2 valves)	50		

Driven by its commitment to advance radiology today and tomorrow, Guerbet Diagnostic Imaging has designed a portfolio of interconnected contrast imaging solutions to enhance your decision-making at each point of the patient journey from diagnosis, to treatment, to follow-up, so you can focus on what matters most, efficiently improving patient outcomes.



UNIK

Tailored interconnected solutions
driving your journey to excellence

www.guerbet.com

Optivantage is a medical device intended for use by qualified healthcare professionals, it is a contrast delivery system and is designed to inject radiopaque contrast media into a patient's vascular system to obtain diagnostic images when used with computed tomography (i.e. «CT») equipment. Class IIb/CE TÜV SÜD 0123 Manufacturer: Liebel-Flarsheim Company LLC EC Rep: Guerbet SA (Villeglante, France)

LF Injector Disposables are a medical devices intended for use by qualified healthcare professionals and designed to be used in the LF contrast injection systems. They contain the contrast and provide the delivery path for the contrast to the patient.

The injection systems are designed to inject contrast media into a patient's vascular system to obtain diagnostic images when used with computed tomography (CT), MR, Angiographic or x-ray equipment. Class IIa/CE TÜV SÜD 0123 Manufacturer: Liebel-Flarsheim Company LLC. Manufacturer: Shenzhen Boon Medical Supply Co., Ltd.

For complete information about precautions and optimal usage conditions, we recommend consulting the instruction for use supplied with the device or by your local Guerbet representative(s).
For use only in countries with applicable health authority registrations.



Scan here
to contact us