

Guerbet Diagnostic Imaging

INTERCONNECTED SOLUTIONS

Driven by its commitment to advance radiology today and tomorrow, **Guerbet Diagnostic Imaging** has designed a portfolio of **interconnected contrast imaging solutions** to enhance your decision-making at each point of the patient journey from diagnosis, to treatment, to follow-up, so you can focus on what matters most, efficiently improving patient outcomes.



HIGH QUALITY CONTRAST MEDIA

For enhanced clinical
decision-making



- Market-leading contrast media for CT, Cath Lab and MRI imaging
- Strong medical support and educational programs

SMART INJECTION SOLUTIONS

For safe and efficient
injections



- Large, flexible range of injectors and disposables
- Tailored to your practice and to patients' needs
 - Reliable technical support services


INNOVATIVE DIGITAL SOLUTIONS

For heightened patient
safety and productivity



- Software Dose&Care®: radiation dose monitoring and tracking
- Software Contrast&Care®: patient injection data and contrast media monitoring and tracking

THE INTERCONNECTED SOLUTIONS ARE NOW COMING TO LIFE FOR YOU



UNIK

Tailored interconnected solutions
driving your journey to excellence

UNIK is a flexible solution, tailored to your specific needs and challenges, designed to add value at each step of the patient healthcare journey, interconnecting our comprehensive portfolio of products, solutions and services in diagnostic imaging.

Contact your Guerbet representative to design the solution
best adapted to your needs and challenges.

www.guerbet.com



Guerbet | 

Dose&Care®

Dose&Care® is a software application which can be integrated into the imaging centre workflow. Dose&Care® records and presents X-ray exposure data and provides a set of tools to document and analyse this information. Dose&Care® is for use by medical imaging and diagnostic health professionals, in particular radiologists, radiographers, medical physicists and personnel responsible for X-ray security and imaging centre administration. Dose&Care® can connect to the centre information system (including PACS, RIS, HIS) or imaging devices and accessories (scanners, X-ray tables, interventional tables, etc.) using DICOM and HL7 communication standards. Dose&Care® collates data into a centralised repository and presents X-ray dose history for each patient, with options to display the information graphically and numerically, and in the form of effective dose (following ICRP 60 and ICRP 103 publications), DLP, CTDI, CTDIvol and SSDE. For comprehensive information and optimal utilisation of Dose&Care, please refer to the User's Manual. Dose&Care® is class I/CE marked Medical device. Legal Manufacturer: MPTronic, 53 rue du Temple 75004 Paris, France. Distributor: Guerbet, BP 57400, 95943 Roissy CdG cedex, France, tel: +33 (0)1 45 91 50 00

Contrast&Care®

Contrast&Care® enables imaging centers to collect, archive, review, and share patient injection data, including contrast media data, adverse events, injector activity, eGFR information and other preexam alerts, such as previously reported allergies. Contrast&Care® also provides options to review protocols, create protocol libraries, and visualize analytics and trends relating to injection activity and contrast media usage. For a comprehensive description of Contrast&Care® and its capabilities, please refer to the user manual. HIPPA and GDPR are taking into account in the design of Contrast&Care®. Contrast&Care® is class I/CE marked Medical device. Legal Manufacturer: MEDEX, 240 allée Jacques Monod, 69800 Saint-Priest. Distributor: Guerbet, BP 57400, 95943 Roissy CdG cedex, France, tel: +33 (0)1 45 91 50 00

For more information, please contact your local Guerbet representative.