

Illumena[®] Néo Single-use
CathLab injector & its
certified* consumables
work better together



Choose the ideal
COMBINATION

* Certified by Guerbet

Essential certified by Guerbet

Combining your **Illumena® Néo** with **Guerbet consumables ensures optimum performance.**

Our priority is to support you by improving our consumables portfolio, thanks to a partnership with selected manufacturers to ensure high quality standards in your daily injection workflow.



All consumables are selected, tested and certified to meet the highest Guerbet quality standards. Legal manufacturers are ISO 13485** accredited and all consumables are CE marked and comply with the European Union regulatory requirements for Medical Devices*** or with the U.S. regulations or with local regulatory requirements when products are approved or registered.

* Certified by Guerbet as compatible on Guerbet injectors - data on file. ** Legal manufacturers Liebel-Flarsheim Company LLC and Shenzhen Boon Medical Supply Co., Ltd. are ISO 13485 certified. *** The EU Medical Device Directive (93/42/EEC, MDD) is a legislation which addresses essential principles for the safety and performance of medical devices during their lifecycle. The new EU Medical Devices Regulation ((EU) 2017/745, MDR) will repeal the existing MDD and bring significant regulatory changes.

Single Angio Syringe with spike

Single-use – 150mL

Pressure limit: 1200psi

(Not compatible with Adaws)

REF #236303




High pressure braided tube

Single-use – 120cm


Pressure limit: 1200psi

REF #223475



Reference number	Designation article	Quantity in box	Requested quantity	Comments
236303	150mL Single Angio Syringe	50		
223475	120cm High pressure braided tube	25		

Driven by its commitment to advance radiology today and tomorrow, Guerbet Diagnostic Imaging has designed a portfolio of interconnected contrast imaging solutions to enhance your decision-making at each point of the patient journey from diagnosis, to treatment, to follow-up, so you can focus on what matters most, efficiently improving patient outcomes.



UNIK
Tailored interconnected solutions
driving your journey to excellence

www.guerbet.com

Ilumena® Néo is a medical device intended for use by qualified healthcare professionals, it is a contrast delivery system and is designed to inject radiopaque contrast media into a patient's vascular system to obtain diagnostic images when used with computed tomography (i.e. «CT») equipment. Class IIb/CE TÜV SÜD 0123 Manufacturer: Liebel-Flarsheim Company LLC. EC Rep: Guerbet SA (Villepinte, France).

LF Injector Disposables are a medical devices intended for use by qualified healthcare professionals and designed to be used in the LF contrast injection systems. They contain the contrast and provide the delivery path for the contrast to the patient.

The injection systems are designed to inject contrast media into a patient's vascular system to obtain diagnostic images when used with computed tomography (CT), MR, Angiographic or x-ray equipment. Class IIa/CE TÜV SÜD 0123 Manufacturer: Liebel-Flarsheim Company LLC. Manufacturer: Shenzhen Boon Medical Supply Co., Ltd.

For complete information about precautions and optimal usage conditions, we recommend consulting the instruction for use supplied with the device or by your local Guerbet representative(s).
For use only in countries with applicable health authority registrations.



**Scan here
to contact us**

Guerbet |